

1 – Final Results - Ladies 500 meter

Pos	Nr	Competitor	Nat	Pair	I/O	Final Time	Points	Behind	Info
1	11	Manli Wang	CHN	17	O	38.31	38.310		
2	18	Jenny Wolf	GER	15	O	38.33	38.330	0.020	
3	39	Svetlana Zhurova	RUS	15	I	38.35	38.350	0.040	
4	29	Sang-Hwa Lee	KOR	12	O	38.37	38.370	0.060	TRC
5	34	Marianne Timmer	NED	14	I	38.59	38.590	0.280	
6	9	Hui Ren	CHN	13	O	38.63	38.630	0.320	
7	24	Sayuri Osuga	JPN	17	I	38.64	38.640	0.330	
	26	Sayuri Yoshii	JPN	18	O	38.64	38.640	0.330	
9	23	Tomomi Okazaki	JPN	16	I	38.70	38.700	0.390	
10	22	Chiara Simionato	ITA	19	I	38.84	38.840	0.530	
11	10	Beixing Wang	CHN	19	O	38.89	38.890	0.580	
12	43	Jennifer Rodriguez	USA	18	I	38.96	38.960	0.650	
13	16	Judith Hesse	GER	11	O	38.99	38.990	0.680	
14	33	Sanne vd Star	NED	2	O	39.04	39.040	0.730	
15	31	Annette Gerritsen	NED	12	I	39.08	39.080	0.770	
16	17	Sabine Völker	GER	14	O	39.18	39.180	0.870	
17	5	Shannon Rempel	CAN	16	O	39.30	39.300	0.990	
18	25	Aki Tonoike	JPN	6	O	39.33	39.330	1.020	
19	27	Seung-Yong Choi	KOR	13	I	39.39	39.390	1.080	
20	35	Marieke Wijsman	NED	7	O	39.46	39.460	1.150	
21	38	Yulia Nemaya	RUS	9	O	39.59	39.590	1.280	
22	44	Amy Sannes	USA	7	I	39.67	39.670	1.360	
23	2	Anzhelika Kotyuga	BLR	11	I	39.68	39.680	1.370	
24	6	Kerry Simpson	CAN	8	O	39.79	39.790	1.480	
25	42	Elli Ochowicz	USA	2	I	39.81	39.810	1.500	
26	28	You-Lim Kim	KOR	1	I	40.06	40.060	1.750	TRC
27	45	Chris Witty	USA	6	I	40.14	40.140	1.830	
28	4	Krisy Myers	CAN	1	O	40.36	40.360	2.050	
29	40	Paulina Wallin	SWE	4	O	40.39	40.390	2.080	
30	8	Kim Weger	CAN	8	I	40.61	40.610	2.300	
31	3	Svetlana Radkevich	BLR	10	O	40.94	40.940	2.630	
32	21	Ágota Toth	HUN	9	I	41.59	41.590	3.280	
33	20	Orsolya Fecskés	HUN	5	O	41.70	41.700	3.390	
34	36	Daniela Dumitru	ROM	3	I	42.00	42.000	3.690	
	41	Yelena Myagkikh	UKR	5	I	42.00	42.000	3.690	
36	37	Daniela Oltean	ROM	4	I	42.01	42.010	3.700	
37	1	Anna Badayeva	BLR	10	I	42.69	42.690	4.380	
	19	Pamela Zoellner	GER	3	O	NS			NS

Referee: Mr. Roland E. Maillard (SUI) Start: 13.00 uur
 Starter: Mr. Matthijs de Vrijer (NED) End: 13.50 uur

Information on ice and track conditions can be found on the last page of the protocol